

DUPLISOLE 24

A-Silicone Duplicating Material



DUPLISOLE 24

Duplicating material used in dental laboratories.

ADVANTAGES


- Extreme accuracy of detail reproduction
- Good elasticity and tear resistance
- Easy to dose and mix
- Unaltered in time even after weeks
- Compatible with plasters, investments and resins
- Fast Setting

CHARACTERISTICS

- Mixing ratio 1:1
- High fluidity
- Low final hardness
- Resistant to high temperature
- Long dimensional stability
- Light blue color

Piece by Piece
with *pleasure.*

Premium Plus Japan commitment
to the own manufacture products



TECHNICAL DATA

Duplisolet 24

Mixing time	1'
Total working time (23°C)	5'
Setting time (23°C)	20'
Hardness after 24 hours	24 Shore A
Linear dimensional change	< 0,05%

PACKAGING

PJ9990	Duplisolet 24 1+1Kg
PJ9990-S	Duplisolet 24 250+250g

MAIN APPLICATION

Duplication of models in fixed
and removable prostheses.

Available by the authorized dealer