

SOLEBITE S

Extra Hard

Bite Registration addition silicone (VPS)





SOLEBITE S

Extra Hard

**Bite registration addition silicone (VPS)
used in dental clinic and laboratory.**

ADVANTAGES

- Very high final hardness
- Easy to be injected on the teeth
- Low bite resistance from the patient
- Quick polymerization
- Easy to see occlusal contact point
- Long dimensional stability
- Easy to trim

CHARACTERISTICS

- Extra Hard (95 Sh-A)
- Extra Fast Setting
- Biocompatible
- Thixotropic
- Dark blue color

Piece by Piece
with *pleasure.*

Premium Plus Japan commitment
to the own manufacture products



TECHNICAL DATA SoleBite S Extra Hard

| | |
|---------------------------|-------------|
| Mixing time | 20" |
| Total working time (23°C) | 30" |
| Time in mouth | 50" |
| Hardness after 24 hours | >95 Shore A |
| Linear dimensional change | < 0,10% |
| Detail reproduction | <20 μm |

PACKAGING

| | |
|-----------|-------------------------------------|
| PJ8190-Ex | Solebite S Ex Hard 2x50ml + 3mixers |
| PJ8190-S | Solebite S Ex Hard 1x50ml + 1mixer |
| DC-0123 | Mixing tips green 50pcs |
| DEL5025 | Bite Tips Eco 50pcs |
| PJ8010 | Dispenser PJ |

MAIN APPLICATION

Occlusal or vestibular bite registration keys of partial / full arch.

Available by the authorized dealer