

# SOLEBITE S

Bite Registration addition silicone (VPS)



## SOLEBITE S

Bite registration addition silicone (VPS)  
used in dental clinic and laboratory.

### ADVANTAGES

- Easy to be injected on the teeth
- Low bite resistance from the patient
- Quick polymerization
- High final hardness
- Easy to see occlusal contact point
- Long dimensional stability
- Easy to trim

### CHARACTERISTICS

- Hard (90 Sh-A)
- Fast Setting
- Biocompatible
- Thixotropic
- 3 years shelf-life
- Light blue color

Piece by Piece  
with *pleasure.*

Premium Plus Japan commitment  
to the own manufacture products



### TECHNICAL DATA

#### SoleBite S

Total working time (23°C)	20"
Time in mouth	30"
Total setting time (23°C)	50"
Hardness after 24 hours	90 Shore A
Linear dimensional change	< 0,13%
Detail reproduction	<20 μm

### PACKAGING

PJ8110-Ex	Solebite S 50mlx2 + mixers
PJ8110-Ex	Sample Solebite S 50ml
DC-0123	Mixing tips green long 50pcs
DEL5025	Bite Tips Eco 50pcs
PJ8010	Dispenser PJ

### MAIN APPLICATION

Occlusal or vestibular bite registration keys of partial / full arch.

Available by the authorized dealer