

SOLEBITE S

Bite Registration addition silicone (VPS)



SOLEBITE S

Bite registration addition silicone (VPS)
used in dental clinic and laboratory.

ADVANTAGES

- Easy to be injected on the teeth
- Low bite resistance from the patient
- Quick polymerization
- High final hardness
- Easy to see occlusal contact point
- Long dimensional stability
- Easy to trim

CHARACTERISTICS

- Hard (90 Sh-A)
- Extra Fast Setting
- Biocompatible
- Thixotropic
- Light blue color

Piece by Piece
with *pleasure.*

Premium Plus Japan commitment
to the own manufacture products



| TECHNICAL DATA | SoleBite S |
|---------------------------|------------|
| Mixing time | 20" |
| Total working time (23°C) | 30" |
| Time in mouth | 50" |
| Hardness after 24 hours | 90 Shore A |
| Linear dimensional change | < 0,10% |
| Detail reproduction | <20 µm |

| PACKAGING | |
|-----------|----------------------------------|
| PJ8110-Ex | Solebite S 2x50ml + 3mixers |
| PJ8110-S | Solebite S 1x50ml + 1mixer green |
| DC-0123 | Mixing tips green long 50pcs |
| DEL5025 | Bite Tips Eco 50pcs |
| PJ8010 | Dispenser PJ |

MAIN APPLICATION

Occlusal or vestibular bite registration keys of partial / full arch.

Available by the authorized dealer