

SOLEBITE S

Super Fast

Bite Registration addition silicone (VPS)



SOLEBITE S

Super Fast

Bite registration addition silicone (VPS)
used in dental clinic and laboratory.

ADVANTAGES

- Extra fast polymerization
- Low residual of silicone into the mixer (short version)
- Easy to be injected on the teeth
- Low bite resistance from the patient
- High final hardness
- Easy to see occlusal contact point
- Long dimensional stability
- Easy to trim

CHARACTERISTICS

- Extra Fast Setting
- Hard (90 Sh-A)
- Biocompatible
- Thixotropic
- 3 years shelf-life
- Indigo blue color

Piece by Piece
with *pleasure.*

Premium Plus Japan commitment
to the own manufacture products



TECHNICAL DATA SoleBite S Super Fast

Total working time (23°C)	15"
Time in mouth	20"
Total setting time (23°C)	35"
Hardness after 24 hours	90 Shore A
Linear dimensional change	< 0,13%
Detail reproduction	<20 µm

PACKAGING

PJ8180-Ex	Solebite S SF 50mlx2 + mixers
PJ8180-Ex	Sample Solebite S SF 50ml
DC-0127	Mixing tips green short 50pcs
DEL5025	Bite Tips Eco 50pcs
PJ8010	Dispenser PJ

MAIN APPLICATION

Occlusal or vestibular bite registration keys of partial / full arch.

Available by the authorized dealer