

# SOLEBITE S

## Super Fast

Bite Registration addition silicone (VPS)



## SOLEBITE S Super Fast

**Bite registration addition silicone (VPS)  
used in dental clinic and laboratory.**

### ADVANTAGES

- Super fast polymerization
- Low residual of silicone into the mixer (short version)
- Easy to be injected on the teeth
- Low bite resistance from the patient
- High final hardness
- Easy to see occlusal contact point
- Long dimensional stability
- Easy to trim

### CHARACTERISTICS

- Super Fast Setting
- Hard (90 Sh-A)
- Biocompatible
- Thixotropic
- Indigo blue color

Piece by Piece  
with *pleasure.*

Premium Plus Japan commitment  
to the own manufacture products



### TECHNICAL DATA SoleBite S Super Fast

Mixing time	15"
Total working time (23°C)	20"
Time in mouth	35"
Hardness after 24 hours	90 Shore A
Linear dimensional change	< 0,10%
Detail reproduction	<20 µm

### PACKAGING

PJ8180-Ex	Solebite S SF 2x50ml + 3mixers
PJ8180-S	Solebite S SF 1x50ml + 1mixer
DC-0127	Mixing tips green short 50pcs
DEL5025	Bite Tips Eco 50pcs
PJ8010	Dispenser PJ

### MAIN APPLICATION

Occlusal or vestibular bite registration keys of partial / full arch.

Available by the authorized dealer