

## SOLEGEL

Catalyst Gel for C-Silicone used in dental laboratories.

### ADVANTAGES

- Accurate dosing through calibrated tip
- Easy to be mixed

### CHARACTERISTICS

- Gel consistency
- Homogeneous formulation
- Mixed with Sole C-Lab Putty 85
- Red color



### MAIN APPLICATION

Used to polymerize Sole C-Lab Putty 85

### TECHNICAL DATA      Solegel with C-Lab 85

Mixing time	30"
Total working time (23°C)	1'30"
Setting time (23°C)	5'30"
Hardness after 24 hours	85 Shore A
Linear dimensional change	< 0,2%
Elastic recovery	> 99%

### PACKAGING

PJ9950 Solegel 60g