

# DUPLISOLE 24

A-Silicone Duplicating Material



## DUPLISOLE 24

Duplicating material used in dental laboratories.

### ADVANTAGES


- Extreme accuracy of detail reproduction
- Good elasticity and tear resistance
- Easy to dose and mix
- Unaltered in time even after weeks
- Compatible with plasters, investments and resins
- Fast Setting

### CHARACTERISTICS

- Mixing ratio 1:1
- High fluidity
- Low final hardness
- Resistant to high temperature
- Long dimensional stability
- Light blue color

Piece by Piece  
with *pleasure.*

Premium Plus Japan commitment  
to the own manufacture products



### TECHNICAL DATA

#### Duplisole 24

Mixing time	1'
Total working time (23°C)	5'
Setting time (23°C)	20'
Hardness after 24 hours	24 Shore A
Linear dimensional change	< 0,05%

### PACKAGING

PJ9990	Duplisole 24 1+1Kg
PJ9990-S	Duplisole 24 250+250g
PJ9995	Duplisole 24 5+5Kg

### MAIN APPLICATION

Duplication of models in fixed  
and removable prostheses.

Available by the authorized dealer