

METAL TRAYS

Stainless steel
Upper and Lower

XS

S

M

L

XL



PERF-LOCK METAL TRAYS

Stainless steel impression trays.

ADVANTAGES

- Extra rigid - no deformation
- Accurate finishing of the surfaces
- Autoclavable at 134°C
- Mechanical retention to the impression material
- Size indicated by number or abbreviation

CHARACTERISTICS

- Stainless steel
- Thick perforated metal
- Packed into single zip bag
- 5 welding spots on the holder
- XS,S,M,L,XL upper and lower sizes available



MAIN APPLICATION

Upper and lower full arch impressions with elastomers or alginates.

| SIZE | UPPER (MM) | LOWER (MM) |
|------|------------|------------|
| XS 1 | 63x45 | 67x57 |
| S 2 | 63x52 | 70x52 |
| M 3 | 72x53 | 74x56 |
| L 4 | 70x58 | 73x63 |
| XL 5 | 78x65 | 77x60 |
| | | |
| | | |
| | | |

PACKAGING

- FC-0220 PERF-LOCK TRAY UPPER XL
- FC-0222 PERF-LOCK TRAY UPPER L
- FC-0224 PERF-LOCK TRAY UPPER M
- FC-0226 PERF-LOCK TRAY UPPER S
- FC-0219 PERF-LOCK TRAY UPPER XS
- FC-0221 PERF-LOCK TRAY LOWER XL
- FC-0223 PERF-LOCK TRAY LOWER L
- FC-0225 PERF-LOCK TRAY LOWER M
- FC-0227 PERF-LOCK TRAY LOWER S
- FC-0229 PERF-LOCK TRAY LOWER XS
- FC-0228 PERF-LOCK TRAY SET