

METAL TRAYS

Stainless steel
Upper and Lower

XS

S

M

L

XL



PERF-LOCK METAL TRAYS

Stainless steel impression trays

ADVANTAGES

- Extra rigid - no deformation
- Accurate finishing of the surfaces
- Autoclavable at 134°C
- Mechanical retention to the impression material
- Size indicated by number or abbreviation

CHARACTERISTICS

- Stainless steel
- Thick perforated metal
- Packed into single zip bag
- 5 welding spots on the holder
- XS,S,M,L,XL upper and lower sizes available



SIZE	UPPER (MM)	LOWER (MM)
XS 1	65x48	68x50
S 2	68x51	71x55
M 3	70x54	75x55
L 4	75x57	76x57
XL 5	80x63	80x60

PACKAGING

- FC-220 Perf-lock Tray Upper XL 5
- FC-222 Perf-lock Tray Upper L 4
- FC-224 Perf-lock Tray Upper M 3
- FC-226 Perf-lock Tray Upper S 2
- FC-219 Perf-lock Tray Upper XS 1
- FC-221 Perf-lock Tray Lower XL 5
- FC-223 Perf-lock Tray Lower L 4
- FC-225 Perf-lock Tray Lower M 3
- FC-227 Perf-lock Tray Lower S 2
- FC-229 Perf-lock Tray Lower XS 1
- FC-228 Perf-lock Tray Set 10pcs

MAIN APPLICATION

Upper and lower full archimpressions with elastomers or alginates.