

GINGISOLE SOFT

Soft version, medium consistency laboratory vinylpolysiloxane used in dental laboratories.

ADVANTAGES

- Low final hardness (40 Shore A)
- High elastic properties and detail reproduction
- It can be finished with the cutter
- For direct and indirect techniques

CHARACTERISTICS

- Mixing ratio 1:1
- Elastic
- Long dimensional stability
- Light pink color



MAIN APPLICATION

Reproduction of the alveolar mucosa on models with implants.

TECHNICAL DATA

Gingisole soft

Mixing time	automix
Total working time (23°C)	2'
Setting time (23°C)	6'
Hardness after 24 hours	40 Shore A
Linear dimensional change	< 0,2%
Elastic recovery	> 99,8%

PACKAGING

PJ9810 Gingisole Soft 2X50ML
+ 6 Mixing tips and oral tips
+ Separator 10ML