

PPJ MIRRORS

Advanced mirrors for dental clinics



PPJ MIRRORS KIT

Advanced mirrors intro kit used in dental clinics

ADVANTAGES

- 3 different advanced mirrors
PPJ1 MIRROR with endo gauge scale
PPJ2 MIRROR single side with vacuum
PPJ2D MIRROR double side with vacuum
- 99% high reflectivity
- High durability of the mirror (30 layers of ceramic)
- Scratch resistance hardness 9H
- Lightweight
- Autoclavable at 134°C

CHARACTERISTICS

- Intro Kit contents:
1pc of PPJ1 MIRROR size 4 single side
1pc of PPJ2 MIRROR size 4 single side with vacuum
1pc of PPJ2D MIRROR size 4 double side with vacuum
1pc of rotatable adapter
1pc of holder diam. 16mm
1pc of holder adapter 11mm
1.8m of silicone tube



TECHNICAL DATA

Size

4

PACKAGING

- PJ6030 PPJ Mirrors Intro Kit gray

MAIN APPLICATION

Providing a reflected image inspecting the teeth and gums

PPJ MIRROR 1


Advanced mirrors used in dental clinics

ADVANTAGES

- Functional endo gauge scale
- 99% high reflectivity
- High durability of the mirror (30 layers of ceramic)
- Scratch resistance hardness 9H
- Lightweight
- Autoclavable at 134°C

CHARACTERISTICS

- 13g weight
- Single side mirror



TECHNICAL DATA

Size

4

PACKAGING

• PJ6000 PPJ Mirror 1 single 3pcs

MAIN APPLICATION

Providing a reflected image inspecting the teeth and gums

PPJ MIRROR 2


Advanced vacuum mirrors used in dental clinics

ADVANTAGES

- Vacuum function
- 99% high reflectivity
- High durability of the mirror (30 layers of ceramic)
- Scratch resistance hardness 9H
- Autoclavable at 134°C

CHARACTERISTICS

- Vacuum connection
- Single side mirror



TECHNICAL DATA

Size

4

PACKAGING

• PJ6060 PPJ Vacuum Mirror 2 single gray 3pcs

MAIN APPLICATION

Providing a reflected image inspecting the teeth and gums

PPJ MIRROR 2D

Advanced vacuum mirrors used in dental clinics

ADVANTAGES

- Vacuum function
- 99% high reflectivity
- High durability of the mirror (30 layers of ceramic)
- Scratch resistance hardness 9H
- Autoclavable at 134°C

CHARACTERISTICS

- Vacuum connection
- Double side mirror



TECHNICAL DATA

Size

4

PACKAGING

• PJ6090 PPJ Vacuum Mirror 2 double gray 3pcs

MAIN APPLICATION

Providing a reflected image inspecting the teeth and gums

PPJ MIRRORS ACCESSORIES

Advanced vacuum mirrors used in dental clinics

DESCRIPTION

Rotatable adapter, mirror holder, 11mm adapter and tube



PACKAGING

- PJ6100 PPJ Mirror 2 Holder + Adapter
- PJ6110 PPJ Mirror 2 Silicone Tube Gray

MAIN APPLICATION

Providing a reflected image inspecting the teeth and gums

Piece by Piece
with *pleasure.*

Premium Plus Japan Co., Ltd.
Japan headquarters
5th floor, 1-1-31, Nishimidorigaoka,
Toyonakashi, Osaka 560-0005, Japan
Distribution center
1-13-5, Haradaminami, Toyonakashi,
Osaka 561-0805, Japan
Tel. +81 6 6845 0066
Fax. +81 6 6845 0064
www.premiumplus.jp
info.japan@premiumplus.jp