

# PRIMSOLE

Addition silicone (VPS) impression material



# PRIMSOLE PUTTY SOFT FAST

Addition silicone (VPS) impression material used in dental clinics.

## ADVANTAGES

- Resistant: high tear strength
- Accurate detail reproduction
- Long dimensional stability
- Excellent elastic recovery
- Safe
- 3-year shelf life

## CHARACTERISTICS

- Soft consistency
- Putty type 0 - ISO 4823
- Quick time in mouth
- Final hardness 65 Shore A
- Biocompatible
- Salmon pink color



## TECHNICAL DATA

|                           |            |
|---------------------------|------------|
| Mixing time               | 30"        |
| Total working time (23°C) | 1'30"      |
| Time in mouth             | 2'30"      |
| Total setting time        | 4'         |
| Linear dimensional change | < 0,2%     |
| Hardness after 24 hours   | 65 Shore A |

## PACKAGING

- PJ8081 PRIMSOLE Putty Soft Fast Set  
250ml + 250ml jars

## MAIN APPLICATION

Master impression taking with both, single and double stage techniques



# PRIMSOLE PUTTY SOFT NORMAL

Addition silicone (VPS) impression material used in dental clinics.

## ADVANTAGES

- Resistant: high tear strength
- Accurate detail reproduction
- Long dimensional stability
- Excellent elastic recovery
- Safe
- 3-year shelf life

## CHARACTERISTICS

- Soft consistency
- Putty type 0 – ISO 4823
- Long time in mouth
- Final hardness 65 Shore A
- Biocompatible
- Orange color



## MAIN APPLICATION

Master impression taking with both, single and double stage techniques

## TECHNICAL DATA

|                           |            |
|---------------------------|------------|
| Mixing time               | 30"        |
| Total working time (23°C) | 2'         |
| Time in mouth             | 3'30"      |
| Total setting time        | 5'30"      |
| Linear dimensional change | < 0,2%     |
| Hardness after 24 hours   | 65 Shore A |

## PACKAGING

- PJ8080 PRIMSOLE Putty Soft Normal Set  
250ml + 250ml jars

# PRIMSOLE LIGHT FAST

Addition silicone (VPS) impression material used in dental clinics.

## ADVANTAGES

- Resistant: high tear strength
- High detail reproduction
- Long dimensional stability
- Excellent elastic recovery
- Safe
- 3-year shelf life
- Combined with advanced mixing tips – easy to lock

## CHARACTERISTICS

- Thixotropic consistency
- Light bodied type 3 – ISO 4823
- Hydrocompatible
- Quick time in mouth
- Biocompatible
- Scannable
- Dark blue color



## TECHNICAL DATA

|                           |            |
|---------------------------|------------|
| Mixing time               | Automixing |
| Total working time (23°C) | 1'30"      |
| Time in mouth             | 2'30"      |
| Total setting time        | 4"         |
| Linear dimensional change | < 0,2%     |
| Detail reproduction       | < 20 µm    |

## PACKAGING

- PJ8031 PRIMSOLE Light Body Fast Set  
2x50ml cartridge + 12 mixing tips gray

## MAIN APPLICATION

Master impression taking with both, single and double stage techniques

# PRIMSOLE LIGHT NORMAL

Addition silicone (VPS) impression material used in dental clinics.

## ADVANTAGES

- Resistant: high tear strength
- High detail reproduction
- Long dimensional stability
- Excellent elastic recovery
- Safe
- 3-year shelf life
- Combined with advanced mixing tips – easy to lock

## CHARACTERISTICS

- Thixotropic consistency
- Light bodied type 3 – ISO 4823
- Hydrocompatible
- Long time in mouth
- Biocompatible
- Scannable
- Light blue color



## TECHNICAL DATA

|                           |            |
|---------------------------|------------|
| Mixing time               | Automixing |
| Total working time (23°C) | 2"         |
| Time in mouth             | 3'30"      |
| Total setting time        | 5'30"      |
| Linear dimensional change | < 0,2%     |
| Detail reproduction       | < 20 µm    |

## PACKAGING

- PJ8030 PRIMSOLE Light Body Normal Set 2x50ml cartridge + 12 mixing tips gray

## MAIN APPLICATION

Master impression taking with both, single and double stage techniques



# PRIMSOLE MEDIUM NORMAL

Addition silicone (VPS) impression material used in dental clinics.

## ADVANTAGES

- Suitable for both monophasic and putty-medium techniques
- Resistant: high tear strength
- High detail reproduction
- Long dimensional stability
- Excellent elastic recovery
- Safe
- 3-year shelf life
- Combined with advanced mixing tips – easy to lock

## CHARACTERISTICS

- Thixotropic consistency
- Medium bodied type 2 – ISO 4823
- Hydrocompatible
- Long time in mouth
- Biocompatible
- Scannable
- Magenta color



## TECHNICAL DATA


|                           |            |
|---------------------------|------------|
| Mixing time               | Automixing |
| Total working time (23°C) | 2"         |
| Time in mouth             | 3'30"      |
| Total setting time        | 5'30"      |
| Linear dimensional change | < 0,2%     |
| Detail reproduction       | < 20 µm    |

## PACKAGING

- PJ8020 PRIMSOLE Medium Body Normal Set 2x50ml cartridge + 12 mixing tips gray

## MAIN APPLICATION

Master impression taking with both, single and double stage techniques



Piece by Piece  
with *pleasure.*

Premium Plus Japan Co., Ltd.  
Japan headquarters  
5th floor, 1-1-31, Nishimidorigaoka,  
Toyonakashi, Osaka 560-0005, Japan  
Distribution center  
1-13-5, Haradaminami, Toyonakashi,  
Osaka 561-0805, Japan  
Tel. +81 6 6845 0066  
Fax. +81 6 6845 0064  
[www.premiumplus.jp](http://www.premiumplus.jp)  
[info.japan@premiumplus.jp](mailto:info.japan@premiumplus.jp)